



## Spare Parts

Alfa Laval Toftejorg™ SaniMidget SB, SaniMagnum SB & SaniMega SB



100000053-EN5 2022-12

Original manual




## Table of contents

---

*The information herein is correct at the time of issue but may be subject to change without prior notice*

1. SaniMidget SB, SaniMagnum SB & SaniMega SB Standard and UltraPure .....	4
---	---



SaniMidget SB, SaniMagnum SB & SaniMega SB  
Standard and UltraPure

Pos	Qty Denomination	Position 1 Connector AISI 316L	Position 2 Rotor PEEK 450G	UltraPure Position 2 Rotor PEEK MG	Position 3 Stator AISI 316L	Position 4 Clip AISI 316L
<b>SaniMidget SB, SaniMagnum SB &amp; SaniMega SB</b>						
-	SaniMidget SB 360° 1" BPE US, TE16B102-x0 .....	TE16B510	TE16B543 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
-	SaniMidget SB 270° UP 1" BPE US, TE16B132-x0 .....	TE16B510	TE16B540 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 360° 1" BPE US, TE16B152-x0 .....	TE16B515	TE16B552 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMidget SB 270° UP 1½" BPE US, TE16B182-x0 .....	TE16B515	TE16B550 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMidget SB 360° 1" ISO 2037, TE16B202-x0 .....	TE16B511	TE16B543 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 360° 1" BPE US, TE16B203-x0 .....	TE16B512	TE16B543 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 270° UP 1" ISO 2037, TE16B232-x0 .....	TE16B511	TE16B540 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 270° UP 1" BPE US, TE16B233-x0 .....	TE16B512	TE16B540 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMagnum SB 360° 1½" ISO 2037, TE17B102-x0 .....	TE17B500	TE17B543 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMagnum SB 270° UP 1½" ISO 2037, TE17B132-x0 .....	TE17B500	TE17B540 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMagnum SB 360° 1½" BPE US, TE17B203-x0 .....	TE17B502	TE17B543 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMagnum SB 270° UP 1½" BPE US, TE17B233-x0 .....	TE16B502	TE17B540 <sup>1)</sup>	TE16B531	TE16B563 <sup>3)</sup>	
	SaniMega SB 360° Clip on 2" BPE US, TE18B102-x0 .....	TE18B500	TE18B543 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 270° UP Clip on 2" BPE US, TE18B132-x0 .....	TE18B500	TE18B540 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 360° HF Clip on 2" BPE US, TE18B142-x0 .....	TE18B500	TE18B544 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 270° UP HF Clip on 2" BPE US, TE18B152-x0 .....	TE18B500	TE18B545 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 360° Weld on 2" BPE US, TE18B203-x0 .....	TE18B502	TE18B543 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 270° UP Weld on 2" BPE US, TE18B233-x0 .....	TE18B502	TE18B540 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 360° HF Weld on 2" BPE US, TE18B263-x0 .....	TE18B502	TE18B544 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
	SaniMega SB 270° UP HF Weld on 2" BPE US, TE18B293-x0 .....	TE18B502	TE18B545 <sup>1)</sup>	TE18B521	TE18B562 <sup>4)</sup>	
<b>SaniMidget SB, SaniMagnum SB &amp; SaniMega SB UltraPure</b>						
	SaniMidget SB 270° UP Clip on 1" BPE US, TE16B132x5 .....	TE16B510	TE16B542 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 360° Clip on 1" BPE US, TE16B102x5 .....	TE16B510	TE16B546 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 270° UP Weld on 1" ISO 2037, TE16B232x5 .....	TE16B511	TE16B542 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 360° Weld on 1" ISO 2037, TE16B202x5 .....	TE16B511	TE16B546 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 270° UP Weld-on 1" BPE US, TE16B233-x5 .....	TE16B512	TE16B542 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	
	SaniMidget SB 360° Weld on 1" BPE US, TE16B203x5 .....	TE16B512	TE16B546 <sup>1)</sup>	TE16B521	TE16B562 <sup>2)</sup>	

<sup>1)</sup> Wear part

<sup>2)</sup> Q-Doc version: 8010005852

<sup>3)</sup> Q-Doc version: 8010005853

<sup>4)</sup> Q-Doc version: 8010005854

**How to contact Alfa Laval**

Contact details for all countries are continually updated on our website.  
Please visit [www.alfalaval.com](http://www.alfalaval.com) to access the information directly.

© Alfa Laval Corporate AB

This document and its contents is owned by Alfa Laval Corporate AB and protected by laws governing intellectual property and thereto related rights. It is the responsibility of the user of this document to comply with all applicable intellectual property laws. Without limiting any rights related to this document, no part of this document may be copied, reproduced or transmitted in any form or by any means (electronic, mechanical, photocopying, recording, or otherwise), or for any purpose, without the expressed permission of Alfa Laval Corporate AB. Alfa Laval Corporate AB will enforce its rights related to this document to the fullest extent of the law, including the seeking of criminal prosecution.